

# Instruction Sheet for UniJet<sup>®</sup> In-Line Filter Kit (MF-8952)

revised 12/17/03

This kit can be used with any BASi Microbore UniJet Column Kit. Its purpose is to provide a disposable, in-line filter between the sample injector and the column inlet.

---

## Inventory

---

Quantity	Description	Part No.
1	5 cm long x 0.12 mm i.d. 316 SS connecting tube	MR-4006
1	PEEK union and two fittings	MR-4403
10	0.5 µm pore size, 1.5 mm diameter frit in plastic ring	MF-8955
1	Instruction sheet	A-1367

---

## Installation

---

1. If your UniJet microbore column is in an operating chromatograph, turn off the pump and wait for the pressure to decrease to zero before breaking any connection.
2. Assemble the parts as shown in F1.




Figure 1.

---

## Filter Replacement

---

1. Turn off the pump and wait until the pressure decreases to zero before breaking any connection.
2. Disconnect the UniJet column at the inlet union.
3. Turn on the pump. The frit should be pushed out by the mobile phase.

UniJet is a trademark of Bioanalytical Systems, Inc.

**BAS**  
Bioanalytical  
Systems, Inc

2701 Kent Avenue  
West Lafayette  
Indiana 47906  
765.463.4527  
Fax 765.497.1102