




CV NOTES



- Bridging Chlorine
- ◐ Chlorine or Ligand
- Nb

SAMPLE: $\text{Nb}_6\text{Cl}_{14}(\text{PPr}_3^n)_4$

MEDIUM: 0.2 M (t-butyl) $_4\text{N}^+\text{PF}_6^-$
IN CH_2Cl_2

CONC: UNKNOWN

RATE: 200 mV/sec

ETRODE: Pt BEAD

REF: Ag/AgCl

MODEL: CV-1A

METAL ATOM CLUSTER COMPOUNDS OFTEN EXHIBIT AN UNUSUAL NUMBER OF CHEMICALLY AND ELECTROCHEMICALLY REVERSIBLE CYCLIC VOLTAMMETRIC WAVES. IN THE PRESENT CASE THE NEUTRAL CLUSTER MAY BE OXIDIZED TO A CATION AND THEN TO A DICATION, OR REDUCED TO AN ANION. CYCLIC VOLTAMMETRY CAN BE ONE OF THE TRANSITION-METAL CHEMIST'S BEST FRIENDS AS ACKNOWLEDGED BY DOUG KLENDWORTH AND PROF. R.A. WALTON OF PURDUE UNIVERSITY WHO SUBMITTED MATERIAL FOR THIS CV NOTE.



2701 Kent Ave
West Lafayette
Indiana 47906