

Determination of Purines in Microdialysates Using UniJet SepStik Columns

1004

Purpose

Determination of adenosine, guanosine and adenine (F1) in microdialysates.




Figure 1. Structures of adenosine (A), guanosine (B) and adenine (C).

In order to separate and sensitively detect purines in microdialysates, a BASi UniJet SepStik microbore column was used. The 1 mm internal diameter increases the concentrations of the eluting purines up to 21-fold compared to standard LC columns.

Existing Methods

LCEC and LCUV with conventional columns.

Conditions

System: BASi Electrochemical Detector Package with a HPLC pump equipped configured for microbore chromatography.

Column: [UniJet SepStik Kit \(BAS P/N MF-8949\)](#).

The packing was ODS 3 μ m silica in a 100 x 1.0 mm bed.

Mobile Phase: 7 mM NaH_2PO_4 containing 3.5% CH_3OH . Adjust pH to 3.04 after adding CH_3OH .

Flow Rate: 80 $\mu\text{L}/\text{min}$.

Detector: A UV detector equipped with a micro-volume cell.

Wavelength: 260 nm

Range: 0.002 AUFS

Rise Time: 1.0 sec.

Detection Limit: 50 pg injected yielded a S/N of 3.

The injection volume was 5 μL .

Sample preparation

Dialysate was directly injected onto the system.




Figure 2. Detection limit test

A: 500pg of each purine standard

B: Blank

C: 50pg of each purine standard

Peak Identification: 1. hypoxanthine
2. adenine
3. inosine
4. guanosine
5. adenosine




Figure 3. Purines in various rat striatum dialysates.

A: Blank (Ringer's solution)

B, C & D: Various dialysate samples

Peak Identification as in F2.